Kawaguchi Support

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# LEVPROT BIOSCIENCE

THE YEAST COMPANY


## Yeast Expression System Free of Prokaryotic DNA Animal-Free Origin

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**Levprot Bioscience** provides animal-free & prokaryotic dna-free recombinant proteins and enzymes. Our technology delivers pure, animal-free proteins without prokaryotic dna contamination, ensuring safety and reliability.

**Levprot Bioscience's precision fermentation** platform enables highly functional proteins and enzymes with essential post-translational modifications for diverse applications: raw materials for biopharma, for in vitro diagnostics, and for foodtech landscape.

**Levprot Bioscience's** streamlined process scales easily, meeting market demands efficiently and cost-effectively.

# BLOCKING AGENTS FOR *IN VITRO* DIAGNOSTICS

## Recombinant k-casein

**Recombinant k-casein** is a non-bovine version of the k-casein fraction of whole bovine milk casein. Like this common bovine one, it is a useful tool for stabilizing amorphous calcium phosphate, and it acts as an organic adhesive and binder for safety matches. The advantages are that k-casein is soluble and offers more consistency and homogeneity than the bovine option, with an animal-free origin. Moreover, as it has been heterologously expressed in yeast cells, it is also free from *Escherichia coli* DNA, RNases and endotoxins. It is supplied in a glycerol-free buffer solution.

Cat. No.	Description	Size
MT10K-G1KCASHA	Recombinant Bovine k-casein RNase-free expressed in yeast. 10 mg/mL.	10 mg


## Recombinant Bovine Albumin (rBSA)

**rBSA** is a non-animal albumin, equivalent to the common Bovine Serum Albumin (BSA), although heterologously expressed in yeast cells. It prevents non-specific adhesion of enzymes to plastic surfaces, and stabilizes proteins during incubation, with even better efficiency than BSA.

rBSA is not only animal-free, but also free from *Escherichia coli* DNA contaminants, RNases and endotoxins. Moreover, it offers far more consistency and homogeneity than its bovine homologous. It can be supplied in glycerol-free buffer solution.

This product is animal-free, RNase-free, endotoxin-free and it is free from *Escherichia coli* DNA contaminants.


Cat. No.	Description	Size
MT10K-G1RBSAHA	rBSA RNase-free expressed in yeast. 20 mg/mL.	10 mg

**Figure 1:** Three different restriction enzymes are used to digest three different plasmids. Plasmids A and B have two restriction sites and plasmid C has a single restriction site. Before and after digestion are indicated as C (control) and E (enzyme), respectively. M indicates DNA marker.



**Results Figure 1:** Same reaction efficiency was observed when substituting the animal-origin BSA of the commercial buffer with recombinant, non-animal rBSA.

**Analysis of non-animal rBSA DNase activity**

**Figure 2:** Analysis of rBSA DNase activity when compared to a standard, commercially available DNase I. Recombinant, non-animal rBSA shows non-DNase activity (<0.001 % activity in presence of calf thymus DNA at 25 °C) without needing any inactivation extra-step.

## Performance of animal-free rBSA blocking ELISA experiments



**Figure 3:** In this experiment, ELISA results were revealed using HRP and TMB as substrate; the specific reaction is directly correlated to the emitted signal at 450 nm. Performance of commercial animal-origin BSA and recombinant BSA expressed in yeast as blocking agents was compared. Three different pAB concentrations were tested (1, 0.5 and 0.1 µg/mL), against a constant antigen concentration (1 µg/mL).

**Results Figure 3:** The higher blocking capacity of non-animal-derived rBSA minimizes non-specific antigen and antibodies bindings to the plate surface.

Hence, a direct correlation is observed between the decrease in signal emitted by the final reaction and the initial concentration of the pAB: less pAB -> fewer specific reactions -> reduced secondary binding -> decreased signal.

On the contrary, when there are non-specific bindings to the plate surface (the case with animal-origin BSA), the signal decreases less, due to undesired reactions taking place.

**Figure 4:** ELISA results were revealed using HRP and TMB as substrate; specific reaction is directly correlated to 450 nm emitted signal. Performances of commercial animal-origin BSA and rBSA expressed in yeast as blocking agents were compared. Three different pAB concentrations were tested (1, 0.5 and 0.1 µg/mL), against a constant antigen concentration (1 µg/mL). In this case, typical sandwich ELISA was not complete, since first blocking step was not carried out.

**Results Figure 4:** When eliminated the first blocking step of classical ELISA experiments, only with the recombinant, non-animal rBSA the signal decreased accordingly to the pAB initial concentration. Due to its around 10 times higher blocking capacity, first blocking ELISA step can be avoided without efficiency loss when using Levprot's recombinant, non-animal rBSA.